



SEE INSIDE FOR

Precision Products and Services

REPAIR SAVINGS | FREE PROMOTIONS | NEW PRODUCTS | SPARE PARTS

THE BEST SERVICE | THE BEST REPAIRS | THE BEST PRICES

Make it an easy choice; try us today!

CANADA'S #1 SPECIALIZED HANDPIECE/SMALL EQUIPMENT SUPPLIER & AUTHORIZED HANDPIECE REPAIR FACILITY

NEW & IMPROVED!

State-of-the-art products



FREE
NO OBLIGATION
ESTIMATES
on all your repairs



TRADE IN!
GET CREDIT
TOWARDS NEW!



PERIODONTIC, ENDODONTIC AND HYGIENE TREATMENT

NSK ENDO-MATE DT Portable Endo Rotary System



- Easy to use flat control panel
- Torque control from 0.1 to 6.5 Ncm
- Depending on the gear ratio of the head selected
- Lightweight comfortable handpiece
- Dual function on/off switch
- Auto-reverse on, Auto-Stop, Auto-reverse off functions
- Compact and portable desktop unit
- Comes with handpiece

Y141120

LIMITED TIME PROMO

BUY 1 NSK VARIOS 970 or 370 Varios System
GET 3 FREE Additional Tips

- Choose between dry operation, bottled solutions or external water supply
- Touch control mode indicator Endo, Perio, General Mode
- Clear power display for 10 power and irrigation flow setting
- Dual water adjustment controls for bottled or external water
- NSK VARIOS 370 & 970 SYSTEM includes 3 Tips: G4, G6, G8

1-866-729-9229

For Best Price & Current Promo

NSK ENDO-MATE TC2 Portable Cordless Endo Rotary System



- Slim Head and neck assures higher visibility
- Large LCD panel and access-friendly key
- Lightweight, ergonomic design
- Smart Auto Reverse & Alarm Function
- 5 programs for different file systems
- Last program settings can be memorized
- 90min. quick charge with non-contact type charger

Y1001036

KaVo - The MULTiflex-ible Dental Air Scaler



Tip-top power and a wide variety of applications is just a SNAP away with the **KaVo SONICflex**. This dental scaler is certainly the instrument with the broadest range of use and minimally invasive in dentistry procedures, including:

- Dental Hygiene
- Periodontics
- Minimally invasive caries therapy
- Cavity shaping and finishing
- Endodontics

KaVo PROPHYflex Air Abrasion System



- Air Polishing System w/New Improved Design
- Removal of stains & bacterial plaque
- Cleaning of enamel surfaces prior to fissure sealing
- Uses MULTiflex coupler

1000.4672 - Black
1006.9926 - Violet
1006.9928 - Green Apple

NSK VARIOS 970 LED Piezo Ultra Sonic System

The All-in-one Self Contained Unit
Y1001166



\$500 OFF

A large tip selection lets you accomplish many procedures:

- HYGIENE
- NON-INVASIVE RESTORATIVE
- ENDO
- PERIO

NSK VARIOS 370 LED Ultra Sonic System

The Compact and Economical Unit

LED System

#Y1003740

NON-OPTIC System

#Y1003742



ENDO SYSTEMS

WE REPAIR AND SELL OEM PARTS FOR THE ORIGINAL J. MORITA ROOT ZX!!

ENDO SYSTEMS | J. MORITA ROOT ZX® II - APEX LOCATOR WORLD'S BEST SELLING UNIT



Root ZX II®
Endodontic Predictability
Apex Locator

- Continuous, accurate measurement
- No zero-adjustment
- Automatic calibration

TriAuto ZX2
Cordless Endo Dental
Handpiece With ApexLocator



- Apex localization and preparation in a cordless handpiece
 - Unique function for mechanical glide path creation
 - Prevention of file breakage and microcracks
 - Apex localization with millimeter precision
 - Rotatable LCD display for left-handed and right-handed users
 - Small instrument head for direct view of the preparation area
- 24-5902800



**Root ZX II with OTR Endo
Handpiece Rotary System**
Canal Preparation Instrument

- Prepare canals while simultaneously taking measurements
 - Motor flexibility: 8 speed settings
 - 5 Automatic safety functions
 - Fully programmable settings
 - Compact and lightweight handpiece design
- 24-5336392
24-5361257

SLOWSPEEDS

1-866-729-9229
For Best Price & Current Promo

SLOWSPEEDS

NSK FX series

- 6 Month Warranty
- Economical with Quality Performance
 - New One Piece Hygienic Design



**NSK FX25 1:1
Push
button Latch**
C1052002



**NSK FX57 4:1
Reduction Screw
In Prophy**
C1054002



**NSK FX65
1:1 Straight
Nosecone**
H1014002

Buy 4 FX Series Handpieces (Mix & Match) - GET 1 FREE!



**NSK FX204 AIR
MOTOR 4 HOLE**
Shorter, Smaller & Lighter
for unsurpassed balance
M1003002



NEW!

KaVo SmARTmatic Series!

- 1 Year Warranty
- Ease and Agility with the more Slender/ Compact Design
 - Fantastic Performance and Design



**KaVo S10K
Straight
Attachment**
10116731



**KaVo S20K
Slowspeed Push
Button Latch**
10116751



**KaVo S31K
Prophy
Screw-in/Snap-on**
10116761

Buy 3 SMARTmatic Handpieces (Mix & Match) - GET 1 FREE MOTOR!

We offer consumables like Surgical
Tubing and have a wide selection
of tips available



**KaVo 181K
Slowspeed Motor**
05355260

FREE



ZOLAR LASER UNIT

ZOLAR TIPS



Join our list of satisfied customers - try us today!

"The team at Precision Dental offer excellent choices of dental and surgical equipment, and their service is second to none!!! Working with this company has been a great experience and I will continue to use their services in the future."

Dr. Tom wierzicki | Calgary, AB

WWW.PRECISIONDENTAL.CA

*All prices and promos are subject to change without notice.



The Perfect solution for all applications!

HIGHSPEED — ENDO — SLOWSPEED

NEW | Premium ELECTROmatic System

- 4.3" Touch Screen for intuitive and easy operation
- Improve your efficiency and increase your confidence in handling complicated cases
- All-in-one with a Prep and Endo mode
- 5 speed programs in prep mode
- Vivid colour touch display for intuitive set-up
- Endo workflows with each up to 10 files
- Preset Endo file systems
- Dual motor option



Dual motor
10128547

**1-866-729-9229
For Best Pricing**



Single motor
10130311



NEW | NLZ Electric System

Now in a more compact size, with higher safety and flexibility for different installation styles. NSK's all new "NLZ" can help any user of air driven handpieces convert to electric driven handpiece and enjoy the full benefits of this smart system without taking up a lot of space.

NLZ Electric System - No Endo Y1003735
NSK safety functions, 'Contra-Check Function' & 'Contra-SAFE Function'

NLZ E Electric & WAVEONE Endo System Y1003736
Efficient Endodontic Treatment thanks to 'Reciprocating Function' and WAVEONE® GOLD files




CALL FOR PRICING

WaveOne is a registered trademark of Dentsply Sirona Endodontics.

**Simple Plug and Play Hookup on most Delivery units.
BONUS FREE HANDPIECE OFFER—WITH ALL ELECTRICS ABOVE!**

SPRING PROMOTION - HANDPIECE TRADE IN



REDEEM THIS VOUCHER WHEN YOU TRADE IN ANY TWO HANDPIECES TOWARDS NEW PURCHASE

REFER US AND RECEIVE




TOWARDS YOUR NEXT ORDER

OVER \$1799 WHEN THEY USE US
REFERENCE DETAILS REQUIRED



WE ALSO OFFER THESE PRODUCTS



Join our list of satisfied customers - try us today!

"With your help we made a very good decision. I'm very happy with the New Handpieces I got from you! Thanks again for your patience."
Dr. Gurmukh Singh Sahi - Calgary, AB

"I have my electrics in and going this week. They are incredible, it is like trading in a Corolla for a Ferrari. Won't be needing air turbines any longer!"
Dr. R Gardiner - Squamish, BC

THE BEST PRICE & ADVICE

Serving Canada-wide since 1998!

WWW.PRECISIONDENTAL.CA

ELECTRIC PRODUCTS



NSK NLX Nano LED Electric System

Y1001361

The NSK NLX Nano LED has been designed to minimize weight and size to improve the balance of instrument and motor in the hand of the clinician, reducing fatigue during prolonged procedures.

Y1001361

NEW Z Series Electric Family MIX & MATCH **BUY 3 GET 1 FREE!**

QUALITY & PERFORMANCE EVOLVES – THE NEXT PHASE BEGINS!



NSK Z45L Surgical 45 Optic Attachment

C1064001

1ST SURGICAL 45 FOR ELECTRIC MOTORS

Ti-Max Z series handpieces feature smaller head and slimmer neck dimensions in the global market, providing a new level of operational visibility.

Accessing otherwise difficult oral regions is made significantly easier by the miniaturized, ergonomic design.



NSK Z95L 1:5 Speed Increaser Attachment

C1034001



NSK Z25L 1:1 Push Button Latch Attachment

C1038001

1-866-729-9229

For Best Price & Current Promo



KaVo ELECTROmatic Systems

Simple Plug and Play Hookup on most Delivery units.



CALL FOR PRICING

COMFORTbase

- 200XDR handpiece sold separately
- The perfect solution for all Highspeed!
- Constant torque allowing for better control and results
- Perfect balance and lightweight

10128539



COMFORTdrive 200XDR Electric Handpiece

10073570

GET 1 FREE HANDPIECE

Ask for details.



CALL FOR PRICING

PLUS

- The Perfect solution for all High and Slowspeed Applications
- LCD Colour Screen and easy operation
- 3 Speed presets

10128541

MASTER Series



KaVo LUX M25L 1:5

10093630

EXPERT Series



KaVo LUX E25L 1:5

10075550



KaVo LUX E20L 1:1 Latch

10075540



KaVo E20C1:1 Latch—Non Optic

10075541

HANDPIECE CARE & MAINTENANCE

EFFICIENT AND COST EFFECTIVE WAY TO INCREASE YOUR HANDPIECE LONGEVITY!

NSK iCARE AUTOMATIC HANDPIECE LUBRICATOR — C2

Y1002796

- Automatically maintains your handpieces to the highest level of care!
- Very efficient and very inexpensive to operate.
- Cleans, Lubricates and Purges up to (4) Handpieces at one time.
- Fits all brands of Highspeeds, Slowspeeds, Electrics and Surgical Handpieces.



FREE PRODUCT OFFER | HUGE VALUE

Ask for details.

MORITA LUBRINA AUTOMATIC HANDPIECE LUBRICATOR

- Saves time & money:
- Uses only 1/3 the lubricant (compared to manual lubrication)
- Services up to 4 handpieces in only 100 seconds
- Protects your equipment, helps keep it running smoothly
- Built-in chuck maintenance system
- Air suction system eliminates mist for clean operation



SURGICAL PRODUCTS

1-866-729-9229
For Best Price & Current Promo



NSK Ti-MAX X450L Surgical 45 Optic Highspeeds

- More Robust Turbine - 21 Watts of Power
- Designed for sectioning and extraction of impacted teeth and 3rd Molars.

(P1078) NSK
(P1079) KaVo
(P1082001) W&H LED
(P1093) Midwest 5 Hole Optic



NSK Ti-MAX X450 Surgical Non-Optic

(P1088) NSK
(P1092) Midwest 4 Hole



BUY 3
get
1 FREE



KaVo SURGtorque

- 19W of cutting power for faster preps
- Side and rear exhaust
- Unique angled head for increased access

S459 L Surgical 45 Optic
(30004808)




S459 C Surgical 45 Non Optic
(30005063) 4 Hole



CONNECTS TO
MOST COUPLERS!



MIDWEST CONNECTION



SURGICAL PIEZO SYSTEM HIGHEST QUALITY @ LOW-END PRICE!!

SAVE
\$4000

NSK Variosurg Surgical System Optic

- No soft tissue damage
 - Safe use and less post-operations care
 - Superior power cutting capability
 - Eliminates heat generation
- Y1001928

- Extractions
- Ridge Split
- Sinus Lift
- Implant Prep
- Bone Harvesting



NSK Variosurg

LED
LIGHT!



KaVo EXPERTsurg Optic LED Implant Rotary System

- One of the World's lightest and smallest micro-motors
- One touch auto calibration feature
- High performance with 80 Ncm torque
- Intuitive programming and colour screen

10083500

*Handpiece Not Included



KaVo SURGmatic S201 L 20:1 LED Surgical Handpiece

10090470



DON'T
MISS OUT!!!
SAVE \$500

Compatible with all
implant handpieces

POWER, SAFETY & ACCURACY
NO COMPROMISE



NSK Surgic Pro Non-Optic

Surgical rotary system includes 20:1 handpiece.
Y1001950

NSK Surgic Pro LED

Surgical rotary system includes 20:1 handpiece.
200-40,000 RPM
High torque accuracy (5-80 NCM)
Y1001949

NSK SMAX SG20 Push Button Surgical Handpiece (ONLY)**

20:1 Ratio W/External & Internal Irrigation. Compatible With All Surgical Motor Systems.
C1010001



NSK SGA-E2S 1:2 SPEED INCREASER 20 DEGREE SURGICAL NOSECONE**

W/ External Water Irrigation
H265001



J.Morita Needles Available



The NSK Surgic Pro is a compact, stylish and fully integrated system incorporating the unique NSK Advanced Handpiece Calibration (AHC). By individually calibrating each handpiece to the micromotor, AHC absolutely guarantees accurate speed and torque control. The wide range of speed and torque settings allows Clinicians to use Surgic Pro with any implant system brand.

HIGHSPEEDS | MIX & MATCH!



KaVo Smart Series
The perfect choice for affordability and reliability!



KaVo Smart Series Push Button Optic
SMARTorque mini S615L Optic (1008.1643)
SMARTorque S619L Optic (1008.1641)

BUY 3 get 1 FREE

Know the difference & **TRY A DEMO TODAY!**

KaVo Master Series

The Premium "Top of the Line" and made to last!



BUY 3 Kavo M8700L/ M8900 Highspeeds

MASTERtorque mini M8700L Optic (3001.0000)
MASTERtorque M8900L Optic (1008.7800)

GET 1 FREE Highspeed



BUY 3 get 1 FREE

Ask about **Limited Time Offer!**



WORLD'S SMALLEST HEAD AND BODY



NSK S-Max Pico Push Button Optic

Ultra Mini Head & Super Slim Body
P114001 NSK
P1142001 KaVo
P1145003 W&H LED

NSK Highspeeds Connect to Most Couplers



NSK Z Series

Push Button Optic Highspeed

Delivers 26 Watts of torque
Light weight and extremely durable.
Titanium Dura Grip that is easier to hold even when wet. Compatible with most couplers.
Interchange with your existing brand.

NSK Style:

- Z800L Mini (P1110)
- Z900L (P1111)
- Z800KL Mini KaVo (P1112)
- Z900KL Standard KaVo (P1113)
- Z800WL Mini W&H LED (P1120)
- Z900WL Standard W&H LED (P1121)

MIX & MATCH

BUY 3 GET 1 FREE!



NSK Ti-Max X Series Push Button Optic

NSK Style:

- X600L (P496)
- KaVo Style:**
X600KL (P858)
- W&H Style:**
X600WHL for LW SERIES (P1094)
X600WLED (P1101001)



Twinpower Push Button Optic



Ultra E Mini Head:
Morita (16-5356970)
KaVo (16-5357497)
NSK (16-5357500)

Standard Head:
Morita (16-5340888)
KaVo (16-5340608)
NSK (16-5340616)

Surgical 45 Optic:
Morita (16-5344670)
KaVo (16-5350522)

Basic:
4 Hole Non Optic (16-5360838)

1-866-729-9229
For Best Price & Current Promo

CHECK FOR CURRENT
2 + 1 OFFER!

COUPLERS - BUY NOW & SAVE \$50.00!



KaVo 460 LED 6 Pin
10073201



J. Morita CP4-LED
165354978




NSK PTL-LED 6 PIN
P1001600



NSK KCL-LED 6 PIN
(KaVo Style) P1005600

HAPPY CUSTOMERS :-)

We strive to offer our customers the best service & repair!



Absolutely the best! It's just like you're right here in Calgary. Speedy great service and excellent prices on the new NSK Handpieces.

Barb, Dr. Peter C. Chan Calgary, AB

Thanks again for your prompt service Rob...always a pleasure doing business with you!

Cheryl Christmas Fall River, NS



HIGHSPEED REPAIRS

TURBINE PROBLEMS?

- Low Torque • Bur Chatter • Seized • Bur Slipping
- High Pitch • Stuck/Broken Bur • Vibration

OTHER PROBLEMS?

- Water/Chip Air Spray • Coupler Leaking
- No Light • Flickering/Dull LED

WE OFFER:

- Factory OEM Turbines • Euro Aftermarket Turbines
- Turbine Rebuilds

ELECTRICS & ATTACHMENTS

We only replace the necessary parts but you're guaranteed optimum performance. Replacement of bearings, rotors, o-rings, seals, gears, springs, cartridge replacement and rebuilds are just some of the things we offer. All services include adjustment, clean, lube and test.

IMPORTANT: If not properly maintained, Electric handpieces can cause burns. So it's important to check the handpiece for normal operation before each procedure (see below).

Did you know we offer consumables such as surgical tubing?




SAFETY CHECK FOR YOUR ELECTRIC HANDPIECE

DOES IT TURN SMOOTHLY?



With handpiece **NOT** attached to motor, insert and manually rotate bur between fingers. Bur should rotate smoothly without a gravelly feel—there should be **NO** turning restriction or wobble!

HOW DOES IT SOUND?



With handpiece attached to motor, run handpiece at a variety of speeds, and stop and start a couple times—there should be **NO** grinding, rattling or clicking noises!

IS THERE EXCESSIVE HEAT?



With handpiece attached to motor, run handpiece for approximately 20 seconds at maximum speed. Immediately stop running and **with caution** touch head and neck to check for anything abnormal—there should be **NO** heat felt or hot smell!

HASSLE FREE—NO OBLIGATION ESTIMATES WITH QUICK TURNAROUND—CALL US TODAY!

"Amazing group of people that service the dental profession in a truly outstanding manner. Their products are only surpassed by the people that work in the firm."

Shore Side Dentistry
Dr. Peter Taylor - Oakville, ON



WE SPECIALIZE IN HANDPIECE REPAIR + SMALL EQUIPMENT

Precision Dental is my go to for fast quality handpiece repairs. I am very satisfied with their level of service and professionalism.

Dr. David Calgary, AB

Just wanted to say a huge "Thank You" for being good to your word. Appreciate you getting our surgical handpiece back to us so promptly.

Elaine, Dr. T. Greene Delta, BC

I've been a Precision customer for many years and they have what it takes to keep my office running smoothly. I appreciate that I can get an honest, friendly and professional experience each time.

Dr. Peter Walford Hornby Island, BC



MIDWEST®



4-48 Hour Turnarounds | Free No Obligation Estimate | Highest Quality Repair



LAB UNIT REPAIR



SURGICAL UNIT REPAIR



SLOWSPEEDS



1.866.729.9229

WE TAKE CARE OF YOUR REPAIR



HIGHSPEED HANDPIECES



ELECTRIC HANDPIECES



ENDO SYSTEM REPAIR

GET FREE PICKUP to all offices where applicable.
WWW.PRECISIONDENTAL.CA